



Microgravity and Neutron
Crystallography

Edward H. Snell

NASA Laboratory for
Structural Biology

11500EJ11

Acknowledgements

Russell Judge, Mark Pusey, Craig Kundrot, Mark van der Woerd, Cindy Barnes, Richard Kephart – NASA Laboratory for Structural Biology, USA

Gloria Borgstahl, Jeff Lovelace – University of Toledo Chemistry Department, USA

Henry Bellamy, Stanford SSRL, USA

John Helliwell, Manchester University, UK

Naomi Chayen, Imperial College, UK

Dean Myles, EMBL, Grenoble, France

Summary of Talk

- Microgravity reduces sedimentation and buoyancy driven convection
- Larger volume crystals of reduced mosaicity result
- Neutrons scatter weakly and large crystals and large fluxes are needed
- The Laue method provides larger flux but needs low mosaicity crystals
- Microgravity can grow larger crystals enabling many neutron experiments

A photograph of the Space Shuttle Columbia in orbit above Earth. The shuttle is oriented vertically, with its nose pointing towards the bottom of the frame. The Earth's surface is visible as a blue and white pattern of clouds and oceans. The shuttle's external tank and solid rocket boosters are clearly visible. Large solar panel arrays are extended from the orbiter. The text "Why grow crystals in microgravity?" is overlaid in the center of the image.

Why grow crystals in
microgravity?





On the ground:

As the solution surrounding the crystal becomes depleted of the growing macromolecule the solution starts to rise due to density differences.

A convective growth plume of solution flows over the crystal face impeding growth and the quality of crystal packing.

Schlieren photograph of a growth plume rising from a lysozyme crystal (pH 4.0, 0.1M sodium acetate, 5% NaCl at 15°C.

M.L. Pusey, J. Cryst. Growth, 122, 1-7, 1992



In microgravity:

Buoyancy driven convection effects are greatly reduced:

- A zone of depleted material is formed around the crystal. This is termed the depletion zone.
- Crystal growth is dominated by diffusion transport allowing more ordered growth producing a more perfect, higher quality crystal.

Sedimentation of the crystals to the bottom of the well does not occur:

- Crystals are suspended in nutrient for a longer time hence they grow to a larger volume.

Experiments are small volume, have low mass and have high potential scientific and commercial return:

- Many experiments can be flown at one time and even a low success rate would still result in many successful results

Microgravity produces larger more perfect crystals


A photograph of a diamond crystal, likely a rough-cut stone, showing iridescent colors. The crystal is roughly octahedral in shape and exhibits a spectrum of colors including yellow, pink, green, blue, and purple, which are characteristic of light interference in diamond. The text "Experimental History" is overlaid in white with a black outline on the central part of the crystal.

Experimental History




- First flight April 1981 using a sounding rocket
- Since then to the year 2000, NASA flew 45 flights on the Space Shuttle Orbiter with only 5 missions dedicated to microgravity
- The first 7 of these 45 missions had no temperature control
- 8 missions delivered experiments to the Mir space station which suffered from temperature variation (days above 30°C), extreme vibrations/accelerations and even the occasional fire!
- Currently the first experiments have been flown on the International Space Station

Over 185 different macromolecular samples have flown. Some are frequent flyers, others have flown only a limited number of times



Microgravity and Macromolecular Crystallography: Kundrot, C.E., Judge, R.A., Pusey, M.L., & Snell, E.H. Crystal Growth and Design. Crystal Growth and Design, 1, 87-99 (2001).



Improvements seen from microgravity samples (same reference as previous figure)

Results from previous missions:

Up to the end of year 2000, 181 different macromolecules have flown:

- 81 flew only once
- 36 flew only twice
- 12 have flown over 10 times, all exhibiting improved results
- Lysozyme has flown 47 times

On average there are only 40 trials per macromolecule, allowing for duplication only 20 experimental conditions.

Success rates are:

- 20% improved overall,
- 35% for those that flew more than once
- up to 60% for those that have flown four times or more.

Iteration



Specific examples of volume increase

Lysozyme STS-51F, 72,
89-M8

γ -interferon STS-26

Bovine insulin STS-37, 43,
49

Luciferase STS-42

Canavalin STS-42

2 domain CD4 STS-48

Factor D STS-50

α -interferon STS-50, 52,
68

Proline isomerase STS-50

Brain prolyl isomerase STS-
50

Human insulin STS-60,
STS-95

STMV STS-65

Apocructacyanin C1 STS-65

Catalase STS-73

Photosystem 1, STS-73

Pike Parvalbumin STS-94
(neutron study)

HIV Protease STS-85

Argumenter of liver
regeneration, STS-73

raf kinase, STS-73

Concanavalin B, STS-77

alcohol dehydrogenase,
STS-78

α -amylase, STS-79, M4

acid phospholrase a2,
Chinese satellite

L-alanine dehydrogenase,
STS-73

Thaumatococcus (ISS)


Types of Crystallization Hardware



PCAM trays with 7 experimental cells (top).
Nine trays are housed in one cylinder (center).
Six cylinders fit in a thermal carrier and are housed in an EXPRESS rack on ISS (bottom).



DCAM trays shown in Single-locker Thermal Enclosure System (STES)



Top view of VDAs housed in a tray

See
<http://crystal.nasa.gov>

The Laue Method

- Polychromatic incident wavelength
 - Exploits more of the available flux – used for neutron diffraction where scattering is poor and flux is much lower than an X-ray source.
- Limited number of exposures needed for a complete data set
 - A series of still exposures are taken with large ϕ steps between images.
- Reflections suffer from overlap
 - Many are illuminated at any one time due to the wide sampling of reciprocal space.
- Needs low mosaicity, high quality crystals
 - Reduce reflection overlap
- For neutrons, need large volume
 - $>1\text{mm}^3$ to achieve good signal to noise in the data.

MOSFLM

HALFSLICE

PEAKSTATS


BEAM-ish

A



E

Studies on perfection, volume and diffraction quality

Zoom Window (300%)



X: 374 Y: 822
H: -5 K: 13 L: -9
FWHM: 0.0107 MAX@: 14.111 d(A): 3.2599 I/Sig: 73.31



/18GBdisk1/data/ssrl_finephi/mgin1/mgin1_2_004.img

Resolution Rings Reflections Beam Center Sectors Zoom Window

Res: 2.09,2.68,3.91,7.68

F


3x

A photograph of the International Space Station (ISS) in orbit above Earth. The station's complex structure, including the central truss and multiple solar panel arrays, is clearly visible against the blue and white background of the planet. The text "What is perfection?" is overlaid in the center of the image.


What is perfection?

Perfection is both good long-range and short-range order in the crystal

- Two types of order within the crystal
- Long-range – over many repeating units
 - Effects in real space large, small in reciprocal space. Increase results in reduced mosaicity and increased signal to noise.
- Short-range – molecule to molecule
 - Effects in real space small, global effects seen in reciprocal space. Increase results in increased resolution and reduced B-factors and diffuse scattering.
- The two types of order are not necessarily linked.



Long-range order -
Mosaicity



(a) Well ordered crystal


(b) Large domains are misaligned with respect to each other.

(c) Smaller domains give Fourier truncation


(d) Variation in a domain

(e) A number of effects combine


From Snell, Borgstahl and Bellamy, *Methods in Enzymology, Macromolecular Crystallography Part C*. Edited by Charlie Carter and Robert Sweet – to be published




(a)



(b)



(c)



(d)


Evidence for the domain makeup of crystals.

X-ray topographs

From Snell, Borgstahl and Bellamy, *Methods in Enzymology, Macromolecular Crystallography Part C*. Edited by Charlie Carter and Robert Sweet – to be published

Each topograph is a greatly magnified image of a reflection. In (a) and (b) the crystal is 1.1 mm by 0.9 mm in projection and defined regions are seen at the different reflections of (a) and (b). Some scattering is also seen on the crystal edges, probably due to mounting. In (c) and (d) the crystal is 1.5 mm by 1.1 mm in projection. In this case an array of domains is seen separated by a boundary layer. The different reflections (c) and (d) illustrate a region in the lower right of the crystal coming into the Bragg diffracting condition at the current orientation.

Original experiments using microgravity



Snell et al., Acta Cryst. D51 099-1102 (1995)

Identical reflections from microgravity and ground grown lysozyme.


Eight times increase in signal to noise.

The larger illuminated volume only accounted for a doubling.

Microgravity 0.0023 degrees, ground 0.0130 degrees.

How is mosaicity related to rocking width ?

Three contributions from the X-ray beam



Mosaicity is the reflection rocking width with these contributions deconvoluted.

From the theory to practice

Deconvoluting the measured data

Spectral

$$|\phi_R| = \sqrt{L^2 \zeta^2 \gamma_h^2 + \gamma_v^2} + \frac{L}{d} \cos \theta_{hkl} \left[\eta + \left(\frac{\delta\lambda}{\lambda} \right) \tan \theta_{hkl} \right]$$

Measured rocking width

Horizontal dispersion

Vertical dispersion

True mosaicity

To measure the mosaicity, η , record data in fine slices, 0.001 degree, minimize vertical and horizontal divergence (synchrotron radiation) and monochromate the beam.

See Bellamy, H. D., Snell, E. H., Lovelace, J., Pokross, M. and Borgstahl, G. E. O. "The High Mosaicity Illusion: Revealing the True Physical Characteristics of Macromolecular Crystals" Acta Cryst. D56, 986-995 (2000).



Short-range order -
Resolution

Good Short-range order.

Mutated Glucose Isomerase X-ray data recorded at 100K.

Long-range order was poor with mosaicity of 0.2 degrees.

Data from beamline 9-2, SSRL, ADSC Quantum IV detector, 120 s exposure (dose mode equivalent at start).


0.96 Å resolution


← Enlarged corner of detector

Overloaded data collected in later low resolution fast pass

1.4 Å resolution

2.4 Å resolution



A photograph of the Space Shuttle Columbia in orbit above Earth. The shuttle is oriented vertically, with its nose pointing towards the bottom of the frame. The orbiter is attached to a white External Tank and two white Solid Rocket Boosters. Large solar panel arrays are deployed from the sides of the External Tank. The Earth's surface is visible below, showing blue oceans and white clouds. The text "Our Experiments and Results" is overlaid in the center of the image.

Our Experiments and Results

Experiments


- Recombinant Human Insulin
 - Grown in the Protein Crystallization Facility (PCF) by temperature reduction on the STS-95 Space Shuttle mission (Launched October 29th, 1998).
- Chicken Egg-white Lysozyme
 - Grown by liquid-liquid diffusion on the STS-95 Space Shuttle mission
- Thaumatin
 - Grown by liquid-liquid diffusion in the Enhanced Gaseous Nitrogen dewar (EGN) on two successive missions to the International Space Station (Sep 11th-Oct 24th, 2000 and Feb 7th-Mar 21st 2001).
- Glucose Isomerase – to be launched
 - Grown on the ground by the batch method.

Insulin


1150EJ11



Temperature controlled method. Grown in Commercial Protein Crystallization Facility (CPCF)




Commercial production for pharmaceutical dosing purposes – uniform morphology and size of crystals from microgravity




Images to same scale,
sedimentation onto the
bottom. Clumping of
crystals.

Microgravity →

Many over 2 mm.
Free floating,
unsedimented.



Data Processing with BEAM-ish



Data processing





Table 1 Diffraction Statistics

Sample	Date [¶]	Orthogonal crystal dimensions (mm)	Crystal Volume (mm ³)	Avg. Max. Intensity (counts)	Avg. η [#] (degrees)	No. Refl.	No. data frames
earth-grown insulin crystals [†]							
earth-1	12/98	0.35X0.35X0.32	0.04	859	0.031 (0.017)	170	2000
earth-2	12/98	0.34X0.26X0.13	0.01	880	0.035 (0.015)	20	500
earth-3	12/98	0.40X0.27X0.19	0.02	914	0.017 (0.005)	174	2000
earth-4	7/99	0.43X0.34X0.19	0.03	202	0.038 (0.024)	14	2000
earth-5	7/99	0.39X0.24X0.22	0.02	590	0.013 (0.004)	172	1999
earth-6	7/99	0.39X0.24X0.17	0.02	431	0.023 (0.010)	72	2000
microgravity-grown insulin crystals [§]							
μ g-1	12/98	0.96X0.88X0.37	0.31	18776	0.004 (0.002)	502	2000
μ g-2	12/98	1.20X0.72X0.48	0.42	19528	0.006 (0.005)	241	1000
μ g-3	12/98	0.90X0.88X0.32	0.25	8195	0.004 (0.004)	176	500
μ g-4	7/99	1.29X0.84X0.43	0.47	12846	0.002 (0.001)	491	2000
μ g-5	7/99	1.72X1.31X0.90	2.04	8362	0.004 (0.002)	489	2000
μ g-6	7/99	1.59X1.59X0.50	1.25	7155	0.003 (0.001)	447	2000



Microgravity Blue


Ground Red

← Structural, short range improvement

Resolution

Physical, long range improvement.

The worst microgravity and best ground crystal data are shown.



Resolution

Summary

- Microgravity crystals
 - Grown by thermal methods
 - had consistently larger diffracting volume > 2 mm in each dimension (34 times larger on average)
 - were more physically perfect (7 times lower mosaicity – 0.004° versus 0.026° on average)
 - (not covered in the talk) cryocooled extremely well for structural data collection
 - Showed a huge improvement in signal to noise

A large, faceted, transparent crystal of Lysozyme is shown, exhibiting iridescent colors. The crystal is roughly octagonal in shape and is set against a dark background. The colors are vibrant and multi-colored, including shades of yellow, pink, green, blue, and orange, which are characteristic of light diffraction through a crystal lattice. The word "Lysozyme" is overlaid in white text on the central part of the crystal.

Lysozyme

Other partitioning results – *i.e.* why we did the experiment

Carter et al., 1999 “Lower dimer impurity incorporation may result in higher perfection of HEWL crystals grown in microgravity A case study”, *J. Crystal Growth* 196, 623-637, report:

A K_{eff} of 9 for ground

A K_{eff} of 2 for microgravity

for a lysozyme dimer impurity in crystallization of lysozyme.

Microgravity was seen to preferentially exclude the dimer – it seemed to act as an impurity filter.

Ground was seen to preferentially include the dimer. A significant result which can easily be tested.

Experimental

- Protein – Lysozyme extracted directly from fresh chicken eggs (Judge et al., Biotechnol. Bioeng. 1998, 59, 776-785).
- Dimer impurity – collected from late eluting fractions of lysozyme preparation on a cation exchange column.
- Crystallization:
 - STS-95, 9 day mission.
 - Microgravity – 0.0%, 0.5%, 0.9%, 1.8% and 3.6% impurities.
 - Ground – same conditions.
- Analysis
 - Size and axial ratios measured
 - Crystals washed and analyzed electrophoretically
 - X-ray data collected at Stanford SSRL beamline 1-5


Experimental – X-ray analysis

Crystal	Dimer	Size(mm)	Course Collection		Superfine Collection		Cell Parameters a=b, c (Å)	Date Collected
			Time	Images	Time	Images		
Microgravity	0%	1.12 x 0.96 x 0.72	20 sec	20	2 sec	2000	78.93, 38.02	7/99
	0.5%	0.62 x 0.61 x 0.37	120 sec	20	5 sec	2000	78.51, 37.50	11/98
	0.9%	0.40 x 0.26 x 0.19	60 sec	20	5 sec	1000	78.87, 37.87	11/98
	1.8%	0.32 x 0.28 x 0.13	120 sec	20	5 sec	2000	79.09, 38.00	11/98
	3.6%	0.54 x 0.37 x 0.26	60 sec	20	5 sec	2000	78.89, 38.06	11/98
Ground	0%	0.86 x 0.40 x 0.40	20 sec	20	2 sec	2000	79.09, 38.09	7/99
	3.6%	0.48 x 0.35 x 0.32	15 sec	20	4 sec	2000	78.87, 38.07	7/99


For each crystal two 10° swathes of coarse data with $\Delta\phi = 1^\circ$ were collected 90° apart. Two 1° swathes of superfine f sliced data as 0.001° separated stills were then collected. For the ground 3.6% case swathes were collected 45° apart. The space group for all crystals was $P4_32_12$.

Results – Lysozyme, short-range improvement


Microgravity,
20s exposure
time




Ground, 20s
exposure time




Microgravity,
0.9% impurity,
60s exposure
time



Microgravity,
3.6% impurity,
60s exposure
time




Crystal Size



Microgravity increases volume, impurities decrease it.

Impurity partitioning



Dimer % (w/w)	3.6	1.8	0.9	0.5
Earth K_{eff}	1.0	1.0	0.7	0.3
Microgravity K_{eff}	0.7	0.6	1.2	2.0


- Within the bounds of error there is no difference in partitioning for 3.6, 1.8 and 0.9% impurities.
- Microgravity preferentially incorporates the dimer at 0.5%

X-ray analysis

		Earth		Microgravity				
		0.0 %	3.6%	0.0 %	0.5 %	0.9 %	1.8 %	3.6 %
Earth	0.0 %	-	22 ref 4.9(0.6)	88 ref 9.9(0.9)	56 ref 9.8(0.4)	42 ref 9.8(0.4)	51 ref 10.0(0.5)	95 ref 9.8(1.0)
	3.6 %	22 ref 9.8(0.4)	-	45 ref 9.4(0.3)	28 ref 9.2(0.3)	13 ref 10.2(0.2)	23 ref 9.9(0.3)	30 ref 9.0(0.2)
Micro-gravity	0.0 %	88 ref 2.1(1.9)	45 ref 0.8(1.0)	-	94 ref 4.0(0.7)	70 ref 1.5(2.1)	102 ref 2.4(2.2)	135 ref 3.7(0.5)
	0.5 %	56 ref 2.1(0.8)	28 ref 2.2(0.8)	94 ref 6.1(0.3)	-	40 ref 1.9(1.0)	74 ref 4.5(0.6)	75 ref 5.1(0.2)
	0.9 %	42 ref 7.7(3.6)	13 ref 6.7(1.8)	70 ref 7.5(7.4)	40 ref 7.3(3.3)	-	51 ref 7.6(4.7)	101 ref 19.7(7.2)
	1.8 %	51 ref 19.6(2.7)	23 ref 5.3(3.8)	102 ref 20.3(5.7)	74 ref 20.2(3.9)	51 ref 20.9(4.2)	-	81 ref 7.6(7.5)
	3.6 %	95 ref 6.3(0.4)	30 ref 1.6(1.1)	135 ref 8.5(0.5)	75 ref 8.1(2)	81 ref 8.2(0.3)	101 ref 8.5(0.3)	-


Symmetry related reflections are compared – Microgravity samples increase in mosaicity with increasing impurity reaching a peak at 1.8%. At 3.6% the mosaicity is lower but still higher than 0% and 5%. Earth showed increased mosaicity with impurity, greater than the equivalent microgravity with the exception of the 1.8% microgravity.

X-ray analysis



Intensity plotted against the initial dimer solution concentration for highest resolution bracket greater than 2.01 Å. Other resolution bins display similar trends. The error bars represent 95% confidence intervals

Signal to noise plot for 0% impurity microgravity and ground crystals.



- The best results (signal/noise, mosaicity and volume) come from using highly purified protein in microgravity.
- Microgravity should not be seen as a step to replace good biochemical practice but may be useful in situations where solution impurities are formed during the crystallization process.
- Microgravity seems to be more sensitive to impurities than the ground.


General notes of caution:

Bigger is not necessarily better

The lysozyme dimer may not be typical of a impurity found in general crystallization experiments.

Our results are the opposite of Carter et al.

Impurity partitioning electrophoretic and mosaicity analysis of microgravity and ground-grown crystals. Snell, E.H et al. *Crystal Growth and Design*, 1, 2, 151-158 (2001).



Thaumatin



Experimental details

- 28 μl 87 mg/ml thaumatin from Sigma degassed and placed in tygon tube then frozen
- 112 μl of 1.25 M sodium potassium tartrate containing 100 mM ADA pH 6.5 was added to the tube and also frozen.
- The tube was sealed
- The tube was allowed to thaw over the duration of the microgravity mission
- X-ray data was collected from beamline 7-1 at SSRL
- Two equal 40 swathes were collected for each crystal studied, 4 microgravity and 2 ground.
- Data was collected in dose mode


	Crystal 1 (HOA) microgravity			Crystal 2 (HOA) microgravity			Crystal 3 (HOH) microgravity		
Size	0.97 x 0.57 x 0.46 mm			0.91 x 0.47 x 0.44 mm			1.97 x 1.00 x 1.00 mm		
Cell	a=b=58.56, c=151.58 Å			a=b=58.53, c=151.59 Å			a=b=58.50, c=151.63 Å		
	Slow	Fast	Combind	Slow	Fast	Combined	Slow	Fast	Combined
Res.	40-1.2	40-1.9	40-1.2	40-1.2	40-1.9	40-1.2	40-1.2	40-1.9	40-1.2
R-factor	9.8(60.1)	5.1(48.1)	7.4(60.1)	7.3(57.1)	5.7(13.0)	8.0(59.0)	6.8(39.0)	5.2(10.6)	6.9(39.0)
I/s	14.2(1.3)	10.7(2.2)	13.9(1.3)	17.3(1.5)	17.2(7.2)	17.5(1.4)	17.6(1.4)	23.9(11.2)	23.8(1.4)
Comp.	64.7(46.1)	42.7(35.6)	74.3(47.3)	85.2(74.5)	70.1(91.1)	89.3(78.2)	83.6(49.3)	86.5(98.8)	90.2(49.3)
Unique ref	54106	9270	61964	71094	15176	74394	69528	18758	75042
Mosaicity	0.122	1.633	0.170	0.106	0.113	0.105	0.090	0.097	0.091

	Crystal 1 (HOO)- ground			Crystal 2 (HOO) - ground			Crystal 4 (IKG) microgravity		
Size	0.24 x 0.16 x 0.16 mm			0.34 x 0.18 x 0.18 mm			0.75 x 0.34 x 0.16 mm		
Cell	a=b=58.49, c=151.42 Å			a=b=58.52c=151.42 Å			a=b=58.56c=151.37 Å		
	Slow	Fast	Combind	Slow	Fast	Combined	Slow	Fast	Combined
Res.	40-1.4	40-1.9	20-1.4	40-1.3	40-1.9	40-1.3	40-1.2	40-1.9	40-1.2
R-factor	5.5(53.3)	10.0(85.3)	13.9(55.6)	7.4(62.9)	11.5(61.7)	15.0(62.6)	4.7(39.4)	3.5(8.6)	5.9(39.3)
I/s	16.8(1.2)	8.7(1.1)	18.6(1.3)	17.1(1.1)	11.2(1.8)	17.5(1.1)	17.9(1.4)	18.4(7.2)	17.8(1.3)
Comp.	89.8(62.5)	87.7(86.4)	90.7(67.3)	85.7(37.8)	91.5(98.1)	86.6(37.7)	72.3(24.3)	55.0(76.6)	72.8(26.1)
Unique ref	47499	18987	47897	56326	196843	56926	60201	11916	60585
Mosaicity	0.102	0.199	0.111	0.118	0.317	0.137	0.148	0.201	0.157

Plots of signal/noise against resolution



← Enlarged plot



Glucose Isomerase – to
be flown

Experimental details

- Grown using the batch method
- 20% (w/w) ammonium sulphate, 50 mg/ml glucose isomerase, pH 7.7, 20°C
- 2.5 mm in longest dimension
- Genetically modified
- Small crystals used for X-ray data collection at SSRL
- Large crystals used for neutron data collection at ILL

Good Short-range order.

Mutated Glucose Isomerase X-ray data recorded at 100K.

Long-range order was poor with mosaicity of 0.2 degrees.

Data from beamline 9-2, SSRL, ADSC Quantum IV detector, 120 s exposure (dose mode equivalent at start).


0.96 Å resolution


← Enlarged corner of detector

Overloaded data collected in later low resolution fast pass

1.4 Å resolution

2.4 Å resolution





Glucose Isomerase:

- Largest crystal sample successfully studied to date by neutron diffraction, 43 kDa, $Z=8$.
- Space group $I222$, cell 93.9 99.7 102.9, resolution 2.5Å
- Complete data sets collected from deuterated and non-deuterated crystals at the ILL, Grenoble, France using the Laue polychromatic technique.
- Judge, Snell, van der Woerd and Myles, work in progress.

A photograph of the International Space Station (ISS) in orbit above Earth. The station's complex structure, including the large solar panel arrays, is clearly visible against the blue and white background of the planet. The word "Summary" is overlaid in white text in the center of the image.

Summary

Summary

- Our experiments in microgravity produced consistently larger volume crystals of insulin, lysozyme and thaumatin.
- Other experiments show the same trend
- Our microgravity experiments showed a consistent reduction in mosaicity
- Our microgravity experiments showed a large increase in signal to noise
- Long range order has been dramatically improved by microgravity. We do not know if short range order is being significantly improved.
- Microgravity can make a big impact on neutron diffraction.

Summary of Talk

- Microgravity reduces sedimentation and buoyancy driven convection
- Larger volume crystals with reduced mosaicity result from microgravity growth
- Neutron scatter weakly and large crystals and large fluxes are needed
- The Laue method provides larger flux but needs low mosaicity
- Microgravity can enable many neutron experiments
- **However:**
 - Neutron diffraction should not be justified just because large perfect crystals are available.
 - Ideally, experiments designed to use the potential of seeing hydrogen's, protonation state or using contrast matching should be aimed for.

Acknowledgements

Russell Judge, Mark Pusey, Craig Kundrot, Mark van der Woerd, Cindy Barnes, Richard Kephart – NASA Laboratory for Structural Biology, USA

Gloria Borgstahl, Jeff Lovelace – University of Toledo Chemistry Department, USA

Henry Bellamy, Stanford SSRL, USA

John Helliwell, Manchester University, UK

Naomi Chayen, Imperial College, UK

Dean Myles, EMBL, Grenoble, France

References to this work

Published

- A test of macromolecular crystallization in microgravity: Large, well-ordered insulin crystals. Borgstahl, G.E.O., et al. *Acta Cryst*, D57, 1204-1207 (2001).
- Impurity partitioning electrophoretic and mosaicity analysis of microgravity and ground-grown crystals. Snell, E.H et al. *Crystal Growth and Design*, 1, 2, 151-158 (2001).
- Microgravity and Macromolecular Crystallography: Kundrot, C.E., et al. *Crystal Growth and Design*. *Crystal Growth and Design*, 1, 87-99 (2000).
- The High Mosaicity Illusion: Revealing the true physical quality of macromolecular crystals: Bellamy, H., et al. *Acta Cryst*. D56, 986-995 (2000).
- BEAM-ish: A graphical user interface for the physical characterization of macromolecular crystals: Lovelace, J., et al. *J. Applied Crystallography*, 33, 1187-1188 (2000).
- Synchrotron X-ray reciprocal space mapping, topography and diffraction resolution studies of macromolecular crystal quality: Boggon, T.J., et al. *Acta Cryst*. D56, 868-880 (2000).
- CCD video observation of microgravity crystallization of lysozyme and correlation with accelerometer data: Snell, E.H., et al. *Acta Cryst*. D53, 747-755 (1997)
- CCD video observation of microgravity crystallization: apocrustacyanin C. Chayen, N.E., et al.,. *J. Cryst. Growth* 171, 219-225 (1997).
- Snell, E.H., et al.. Improvements in lysozyme protein crystal perfection through microgravity growth. *Acta Cryst*. D51, 1099-1102 (1995).

In Preparation

- Macromolecular Crystal Quality, Snell, E.H. et al., *Methods in Enzymology, Macromolecular Crystallography Part C* edited by Carter and Sweet – in preparation
- Thaumatin crystallization aboard the International Space Station using liquid-liquid diffusion in the enhanced gaseous nitrogen dewar (EGN). Barnes, C.L. et al., Submitted to *Acta Cryst D*.

Biological Crystal Growth in Space



Technical Interest

Working with NASA

General Interest

Research Announcements


UPDATED
MAY 29, 2001

Contact Us

News

Site Map

MSFC Only
IBC
Project Team



NASA Sites

Microgravity Research Program

Marshall Space Flight Center

NASA

To fly samples go to <http://crystal.nasa.gov> - an application form is in the technical interest section